

IL18R1/IL-18 R alpha Rabbit pAb

货号: AYP12430

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Homo sapiens WB: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	62kDa
实测分子量	75kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NCI-H460
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 19-320 of human IL18R1/IL-18 R alpha (NP_003846.1).
序列	AESCTSRPHITVVEGEPFYLKHCSCSLAHEIETTTKSWYKSSGSQEHVELNPRSSSRIALHDCVLEFWPVELNDTGSYFFQM KNYTQKWKLNVIRRNKHSCFTEQVTSKIVEVKKFFQITCENSYYQTLVNSTSLYKNCKLLENNKNPTIKKNAEFEDQGY YSCVHFLHHNGKLFNITKTFNITIVEDRSNIVPVLGPKLNHVAVELGKNVRLNCSALLNEEDVIYWMFGEENGSDPNIHEEK EMRIMTPEGKWHASKVLRNIENIGESNLNVLYNCTVASTGGTDTKSFILVRKADMA

靶点信息

研究背景	The protein encoded by this gene is a cytokine receptor that belongs to the interleukin 1 receptor family. This receptor specifically binds interleukin 18 (IL18), and is essential for IL18 mediated signal transduction. IFN-alpha and IL12 are reported to induce the expression of this receptor in NK and T cells. This gene along with four other members of the interleukin 1 receptor family, including IL1R2, IL1R1, ILRL2 (IL-1Rrp2), and IL1RL1 (T1/ST2), form a gene cluster on chromosome 2q. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	8809
基因名	IL18R1
Swiss	Q13478
别名	IL18R1;CD218a;CDw218a;IL-1Rrp;IL18RA;IL1RRP

产品验证



Western blot analysis of IL18R1/IL-18 R alpha expressed in NCI-H460 using IL18R1/IL-18 R alpha Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120 S.

实验步骤

访问官网浏览详情: www.ablybio.cn