

— ABLYBIO, Help Your Research



ISG20 Rabbit pAb

货号: **AYP12537**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens WB: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100
理论分子量	9kDa/20kDa
实测分子量	20kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,A-549,MCF7,Mouse spleen,Rat spleen
细胞定位	Cajal body,Cytoplasm,Nucleus,nucleolus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human ISG20 (N P_002192.2).
------	---

靶点信息

研究背景	Interferon-induced antiviral exoribonuclease that acts on single-stranded RNA and also has minor activity towards single-stranded DNA. Exhibits antiviral activity against RNA viruses including hepatitis C virus (HCV), hepatitis A virus (HAV) and yellow fever virus (YFV) in an exonuclease-dependent manner. May also play additional roles in the maturation of snRNAs and rRNAs, and in ribosome biogenesis.
基因ID	3669
基因名	ISG20
Swiss	Q96AZ6 (https://www.uniprot.org/uniprotkb/Q96AZ6/entry)
别名	ISG20,CD25,HEM45,ISG20 Rabbit pAb,Estrogen-regulated transcript 45 protein,Promyelocytic leukemia nuclear body-associated protein ISG20

产品验证



Western blot analysis of ISG20 expressed in Mouse spleen, Rat spleen, Raji using ISG20 Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)