

NUP107 Rabbit pAb

货号: AYP12582

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IF: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	79kDa/103kDa/106kDa
实测分子量	107kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,HeLa,293T,Jurkat,HepG2,Mouse brain,Mouse kidney
细胞定位	Chromosome,Nucleus,Nucleus membrane,centromere,kinetochore,nuclear pore complex
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 766-925 of human NUP 107 (NP_065134.1).
序列	HMNSVPQKPALIPQPTFTEKVAHEHKEKKYEMDFGIWKGHLTDALADVKEKMYNVLLFVDGGWMVDVREDAKEDHERT HQMVLLRKLCLPMLCFLHLHTLHSTGQYQECLQLADMVSSERHKLYLVFSKEELRKLKQLRESSMLLDQGLDPLGYEIQL

靶点信息

研究背景	This gene encodes a member of the nucleoporin family. The protein is localized to the nuclear rim and is an essential component of the nuclear pore complex (NPC). All molecules entering or leaving the nucleus either diffuse through or are actively transported by the NPC. Alternate transcriptional splice variants of t his gene have been observed but have not been thoroughly characterized.
------	--

基因ID	57122
基因名	NUP107
Swiss	P57740
别名	NUP107;NPHS11;NUP84

产品验证



Western blot analysis of NUP107 expressed in Mouse brain, BxPC-3, HeLa using NUP107 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn