

Snail Rabbit pAb

货号: **AYP12614**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	29kDa
实测分子量	29kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,NIH/3T3,Mouse lung,Mouse heart,Rat lung
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human Snail (O 95863).
序列	MPRSFLVRKPSDPNRKPNYSELQDSNPEFTFQQPYDQAHLLAAIPPEILNPTASLPLIWDVSVLAPQAQPIAWASLRLQES PRVAELTSLSDEDSGKGSQPPSPAPSSFSSTSVSSLEAEYAAFPGLGQVPKQLAQLSEAKDLQARKAFNCKYCNKEY LSLGALKMHIRSHLPCVCGTCGKAFSRPWLLQGHVTRHTTGEKPFSCPHCSRAFAADRSNLRAHLQTHSDVKKYQCQACA RTFSRMSLLHKHQESGCSGCPR

靶点信息

研究背景	The Drosophila embryonic protein snail is a zinc finger transcriptional repressor which downregulates the expression of ectodermal genes within the mesoderm. The nuclear protein encoded by this gene is structurally similar to the Drosophila snail protein, and is also thought to be critical for mesoderm formation in the developing embryo. At least two variants of a similar processed pseudogene have been found on chromosome 2.
基因ID	6615
基因名	SNAI1
Swiss	O95863
别名	SNAI1;SLUGH2;SNA;SNAH;SNAIL;SNAIL1;dj710H13.1

产品验证



Western blot analysis of Snail expressed in Mouse lung, Rat lung, HeLa using Snail Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn