

RBBP8 Rabbit pAb

货号: AYP12666

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IF: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	98kDa/101kDa/102kDa
实测分子量	110KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat thymus
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 648-897 of human RBBP 8 (NP_002885.1).
序列	VSFENIQWSIDPGADLSQYKMDVTVIDTKDGSQSKLGGETVDM DCTLVSETVLLKMKKQEQKGEKSSNEERKMNDSLE DMFDRTTHEEYESCLADSFSAADDEEELSTATKKLHTHGDKQDKVKQKAFVEPYFKGDERETSLQNFPHIEVVRKKEER RKLLGHTCKECEIYYADMPAEEREKKLASC SRHRFRYIPPNTPENFWEVGFPSTQTCMERYIKEDLDPCPRPKRRQPYNAI FSPKGKEQKT

靶点信息

研究背景	The protein encoded by this gene is a ubiquitously expressed nuclear protein. It is found among several proteins that bind directly to retinoblastoma protein, which regulates cell proliferation. This protein complexes with transcriptional co-repressor CTBP. It is also associated with BRCA1 and is thought to modulate the functions of BRCA1 in transcriptional regulation, DNA repair, and/or cell cycle checkpoint control. It is suggested that this gene may itself be a tumor suppressor acting in the same pathway as BRCA1. Three transcript variants encoding two different isoforms have been found for this gene. More transcript variants exist, but their full-length natures have not been determined.
基因ID	5932
基因名	RBBP8
Swiss	Q99708
别名	RBBP8;COM1;CTIP;JWDS;RIM;SAE2;SCKL2

产品验证



Western blot analysis of RBBP8 expressed in Rat thymus using RBBP8 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn