

— ABLYBIO, Help Your Research



# PTGES Rabbit pAb

货号: **AYP12774**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>IHC:</b> Homo sapiens <b>WB:</b> Sus scrofa
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	
实测分子量	17kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Mouse lung,Rat lung
细胞定位	Membrane
纯化	Affinity purification


## 抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PTGES (NP_004869.1).
------	--

## 靶点信息

研究背景	
基因ID	9536
基因名	PTGES
Swiss	O14684 ( <a href="https://www.uniprot.org/uniprotkb/O14684/entry">https://www.uniprot.org/uniprotkb/O14684/entry</a> )
别名	PTGES Rabbit pAb,PTGES,Glutathione peroxidase PTGES,Glutathione transferase PTGES,Microsomal glutathione S-transferase 1-like 1,Microsomal prostaglandin E synthase 1,p53-induced gene 12 protein,MGST1 L1,MPGES1,PGES

## 产品验证



Western blot analysis of PTGES expressed in MCF7, Mouse lung, Rat lung using PTGES Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)