

IL33 Rabbit pAb

货号: **AYP12923**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | WB: Human corneal epithelial cells , Mus musculus |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:1000 |
| 理论分子量 | 16kDa/25kDa/26kDa/30kDa |
| 实测分子量 | 33KDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | HepG2,A-549,Rat lung,Mouse lung,Mouse kidney,Mouse stomach |
| 细胞定位 | Chromosome,Cytoplasmic vesicle,Nucleus,Secreted,secretory vesicle |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 110-190 of human IL33 (NP_254274.1). |
| 序列 | DSSITGISPITEYLASLSTYNDQSITFALEDESYEIYVEDLKKDEKKDKVLLSYYESQHPNESGDGVDGKMLMVTLSPTK |

靶点信息

| | |
|------|--|
| 研究背景 | The protein encoded by this gene is a cytokine that binds to the IL1RL1/ST2 receptor. The encoded protein is involved in the maturation of Th2 cells and the activation of mast cells, basophils, eosinophils and natural killer cells. Several transcript variants encoding different isoforms have been found for this gene. |
| 基因ID | 90865 |
| 基因名 | IL33 |

| | |
|--------------|---|
| Swiss | O95760 |
| 别名 | IL33;C9orf26;DVS27;IL1F11;NF-HEV;NFEHEV |

产品验证



Western blot analysis of IL33 expressed in Mouse lung,Rat lung,HepG2 using IL33 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn