

UIMC1 Rabbit pAb

货号: **AYP12925**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | IF: human cells , Mus musculus WB: Mus musculus |
| 应用 | WB IF/ICC |
| 推荐浓度 | WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100 |
| 理论分子量 | 17kDa/38kDa/61kDa/70kDa/79kDa |
| 实测分子量 | 120kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | MCF7,293T,Mouse thymus,Rat brain |
| 细胞定位 | Nucleus |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human UIMC1 (NP_057374.3). |
| 序列 | MPRRKKKVKVESESRNLEKKDVETTSSVSVKRKRRLLEDAFIVISDSGDGEEPKEENGLQKTKTKQSNRAKCLAKRKIAQMTE EEQFALALKMSEQEAREVNSQEEEEELLRKAIAESLNSCRPSDASATRSRPLATGPSSQSHQEKTDS |

靶点信息

| | |
|------|--|
| 研究背景 | This gene encodes a nuclear protein that interacts with Brca1 (breast cancer 1) in a complex to recognize and repair DNA lesions. This protein binds ubiquitinated lysine 63 of histone H2A and H2AX. This protein may also function as a repressor of transcription. Alternative splicing results in multiple transcript variant s. |
|------|--|

| | |
|-------|-----------------------|
| 基因ID | 51720 |
| 基因名 | UIMC1 |
| Swiss | Q96RL1 |
| 别名 | UIMC1;RAP80;X2HRIP110 |

产品验证



Western blot analysis of UIMC1 expressed in Mouse thymus, Rat brain, MCF7 using UIMC1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn