

# Phospho- $\beta$ -Catenin-S33/S37/T41 Rabbit pAb

货号: **AYP12992**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Mus musculus
应用	<a href="#">IHC</a>
推荐浓度	<b>IHC:</b> 1:50 - 1:200
理论分子量	85kDa
实测分子量	85kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HEK-293 cells, NIH/3T3 cells, HeLa cells, mouse brain tissue, A431 cells, MCF-7 cells, mouse liver tissue, rat brain tissue, C6 cells, rat liver tissue
细胞定位	Cell junction,Cell membrane,Cytoplasm,Nucleus,adherens junction,centrosome,cytoskeleton,microtubule organizing center,spindle pole
纯化	Affinity purification

## 抗原信息

抗原信息	A synthetic phosphorylated peptide around S33 & T41 of human beta Catenin (NP_001091679).
序列	SGIHSGATT

## 靶点信息

研究背景	The protein encoded by this gene is part of a complex of proteins that constitute adherens junctions (AJs). AJs are necessary for the creation and maintenance of epithelial cell layers by regulating cell growth and adhesion between cells. The encoded protein also anchors the actin cytoskeleton and may be responsible for transmitting the contact inhibition signal that causes cells to stop dividing once the epithelial sheet is complete. Finally, this protein binds to the product of the APC gene, which is mutated in adenomatous polyposis of the colon. Mutations in this gene are a cause of colorectal cancer (CRC), pilomatixoma (PTR), medulloblastoma (MDB), and ovarian cancer. Alternative splicing results in multiple transcript variants.
------	--

基因ID	1499
基因名	CTNNB1
Swiss	P35222
别名	CTNNB;MRD19;armadillo;beta Catenin;CTNNB1

产品验证

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)