

MALT1 Rabbit pAb

货号: **AYP13027**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	91kDa/92kDa
实测分子量	102kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T
细胞定位	Cytoplasm,Nucleus,perinuclear region
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 125-450 of human MALT 1 (NP_776216.1).
序列	PGIKITVNPESKAVLAGQFVKLCCRATGHPFVQYQWFKMNKEIPNGNTSELIFNAVHVKDGAFYVCRVNNNFTFEFSQWS QLDVCDIPESFQRSVDGVSESKLQICVEPTSQKLMPGSTLVLQCVAVGSPIPHYQWFKNELPLTHETKKLYMVPYVDLEHQ GTYWCHVYNDRDSQDSKKVEIIIDELNNLGHDPNKEQTTDQPLAKDKVALLIGNMNYREHPKLKAPLVDVYELTNLLRQLD FKVWSLLDLTEYEMRNAVDEFLLLLLDKGVYGLLYAGHGYENFGNSFMVPVDAPNPYRSENCLCVQNILKLMQEKETGLNV FLL

靶点信息

研究背景	This gene has been found to be recurrently rearranged in chromosomal translocation with two other genes - baculoviral IAP repeat-containing protein 3 (also known as apoptosis inhibitor 2) and immunoglobulin heavy chain locus - in mucosa-associated lymphoid tissue lymphomas. The protein encoded by this gene may play a role in NF-kappaB activation. Two alternatively spliced transcript variants encoding different isoforms have been described for this gene.
基因ID	10892
基因名	MALT1
Swiss	Q9UDY8
别名	MALT1;IMD12;MLT;MLT1;PCASP1

产品验证



Western blot analysis of MALT1 expressed in HeLa,293T using MALT1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn