

— ABLYBIO, Help Your Research



LASP1 Rabbit pAb

货号: **AYP13038**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IF/ICC IP
推荐浓度	WB: 1:1000 - 1:3000 IF/ICC: 1:50 - 1:200 IP: 1:500 - 1:1000
理论分子量	23kDa/29kDa/36kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,BT-474,SKOV3,U-87MG,LO2,Mouse liver,Mouse brain
细胞定位	Cytoplasm,cell cortex,cytoskeleton
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 130-205 of human LASP 1 (NP_006139.1).
------	--

靶点信息

研究背景	This gene encodes a member of a subfamily of LIM proteins, characterized by a LIM motif and a domain of Src homology region 3, and also a member of the nebulin family of actin-binding proteins. The encoded protein is a cAMP and cGMP dependent signaling protein and binds to the actin cytoskeleton at extensions of the cell membrane. The encoded protein has been linked to metastatic breast cancer, hematopoietic tumors such as B-cell lymphomas, and colorectal cancer.
基因ID	3927
基因名	LASP1
Swiss	Q14847 (https://www.uniprot.org/uniprotkb/Q14847/entry)
别名	LASP1,Lasp-1,MLN50,LASP1 Rabbit pAb,Metastatic lymph node gene 50 protein

产品验证



Western blot analysis of LASP1 expressed in Mouse liver,A-549,BT-474 using LASP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)