

PDCD2L Rabbit pAb

货号: **AYP13062**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens IHC: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	39kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,MCF7,SKOV3,BT-474,Mouse testis,Mouse brain,Mouse heart
细胞定位	cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-358 of human PDCD2 L (NP_115722.1).
序列	MAAVLKPVLLGLRDAPVHGSPTGPGAWTASKLGGIPDALPTVAAPRPVCQRCGQPLALVVQVYCPLEGSPFHRLHVFAC ACPGCSTGGARSWKVFRSQCLQVPEREAQDAQKQGNLSAAEDWCEGADDWGS DTEEGPSPQFTLDFGNDASSAKDV DWTARLQDLRLQDAVLGAAHPVPPGLPLFLPYIICVADEDDYRDFVNL DHAHSLLRDYQQREGIAMDQLLSQSLPNDGD EKYEKTIKSGDQTFYKFMKRIAACQEILRYSWSGEPFLTCTPTSEVTELPAC SQCGGQRIFEFQLMPALVSM LKSANLGLS VEFGTILVYTCEKSCWPPNHQTPMEEFCIIQEDPDELLFK

靶点信息

研究背景	Over-expression suppresses AP1, CREB, NFAT, and NF-κB transcriptional activation, and delays cell cycle progression at S phase.
------	---

基因ID	84306
基因名	PDCD2L
Swiss	Q9BRP1
别名	PDCD2L

产品验证



Western blot analysis of PDCD2L expressed in Mouse testis, HeLa, MCF7 using PDCD2L Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn