

# TNFAIP3 Rabbit pAb

货号: **AYP13122**

## 产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>IF:</b> Mus musculus
应用	<a href="#">WB</a> <a href="#">IF/ICC</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	89kDa
实测分子量	89kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm,Lysosome,Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human TNFAIP 3 (NP_006281.1).
序列	MAEQVLPQALYLSNMRKAVKIRERTPEDIFKPTNGIIHFKTMHRYTLEMFRTCQFCPQFREIIHKALIDRNIQATLESQKCLN WCREVRKLVALKTNGDGNCLMHATSQYMWGVQDLDLVLKALFSTLKETDTRNFKFRWQLESLSQEFVETGLCYDTR NWNDEWDNLIKMASTDTPMARSGLQYNS

## 靶点信息

研究背景	This gene was identified as a gene whose expression is rapidly induced by the tumor necrosis factor (TNF ). The protein encoded by this gene is a zinc finger protein and ubiquitin-editing enzyme, and has been shown to inhibit NF-kappa B activation as well as TNF-mediated apoptosis. The encoded protein, which has both ubiquitin ligase and deubiquitinase activities, is involved in the cytokine-mediated immune and inflammatory responses. Several transcript variants encoding the same protein have been found for this gene .
------	--

基因ID	7128
基因名	TNFAIP3
Swiss	P21580
别名	A20;AISBL;OTUD7C;TNFA1P2;TNFAIP3

## 产品验证



Western blot analysis of TNFAIP3 expressed in HUVEC, rat brain using TNFAIP3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)