

— ABLYBIO, Help Your Research



# SERPINB5 Rabbit pAb

货号: **AYP13146**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Mus musculus
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	25kDa/42kDa
实测分子量	42KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Secreted,extracellular space
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 21-189 of human SERPINB5 (NP_002630.2).
------	---

## 靶点信息

研究背景	Tumor suppressor. It blocks the growth, invasion, and metastatic properties of mammary tumors. As it does not undergo the S (stressed to R (relaxed conformational transition characteristic of active serpins, it exhibits no serine protease inhibitory activity.
基因ID	5268
基因名	SERPINB5
Swiss	P36952 ( <a href="https://www.uniprot.org/uniprotkb/P36952/entry">https://www.uniprot.org/uniprotkb/P36952/entry</a> )
别名	SERPINB5,PI5,maspin,serpin B5,SERPINB5 Rabbit pAb,Peptidase inhibitor 5

## 产品验证



Western blot analysis of SERPINB5 expressed in HeLa using SERPINB5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))