

SFRP5 Rabbit pAb

货号: **AYP13166**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	35kDa
实测分子量	35-38KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	RPMI 8226,A-549,Rat skeletal muscle
细胞定位	extracellular space
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 30-317 of human SFRP5 (NP_003006.2).
序列	EEYDYYGWQAEPLHGSRYSKPPQCLDIPADLPLCHTVGYKRMRLPNLLEHESLAEVKQQAASSWLPLLAKRCHSDTQVFLC SLFAPVCLDRPIYPCRSCLCEAVRAGCAPLMEAYGFPWPEMLHCHKFPLDNDLCIAVQFGHLPATAPPVTKICAQCEMEHSA DGLMEQMCSSDFVVKMRIKEIKIENGDRKLIQAQKKKLLKPGPLKRKDTKRLVLHMKNAGCPCPQLDSLGSFLVMGR KVDGQLLLMAVYRWDKKNKEMKFAVKFMFSYPCSLYYPFFYGAAEPH

靶点信息

研究背景	Secreted frizzled-related protein 5 (SFRP5) is a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. SFRP5 and SFRP1 may be involved in determining the polarity of photoreceptor cells in the retina. SFRP5 is highly expressed in the retinal pigment epithelium, and moderately expressed in the pancreas.
基因ID	6425
基因名	SFRP5
Swiss	Q5T4F7
别名	SFRP5;SARP3

产品验证



Western blot analysis of SFRP5 expressed in RPMI 8226, A-549, Rat skeletal muscle using SFRP5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn