

— ABLYBIO, Help Your Research



SLC30A1 Rabbit pAb

货号: **AYP13173**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Horabagrus nigricollaris
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	55kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human SLC30A1 (NP_067017.2).
------	--

靶点信息

研究背景	May be involved in zinc transport out of the cell.
------	--

基因ID	7779
基因名	SLC30A1
Swiss	Q9Y6M5 (https://www.uniprot.org/uniprotkb/Q9Y6M5/entry)
别名	SLC30A1,ZNT1,ZRC1,SLC30A1 Rabbit pAb,Solute carrier family 30 member 1,Zinc transporter 1

产品验证



Western blot analysis of SLC30A1 expressed in Mouse kidney using SLC30A1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)