



IER3 Rabbit pAb

货号: **AYP13217**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:100
理论分子量	16kDa
实测分子量	17KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,MCF7,Mouse testis
细胞定位	Membrane,Single-pass type II membrane protein
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human IER3 (NP_003888.2).
------	--

靶点信息

研究背景	This gene functions in the protection of cells from Fas- or tumor necrosis factor type alpha-induced apoptosis. Partially degraded and unspliced transcripts are found after virus infection in vitro, but these transcripts are not found in vivo and do not generate a valid protein.
基因ID	8870
基因名	IER3
Swiss	P46695 (https://www.uniprot.org/uniprotkb/P46695/entry)
别名	IER3,DIF-2,DIF2,GLY96,IEX-1,IEX-1L,IEX1,PRG1,IER3 Rabbit pAb,Differentiation-dependent gene 2 protein,Immediate early protein GLY96,Immediate early response 3 protein,PACAP-responsive gene 1 protein

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)