

— ABLYBIO, Help Your Research



# IL20RB Rabbit pAb

货号: **AYP13276**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>neutralizing:</b> Homo sapiens
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	16kDa/35kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,HeLa,HepG2,Mouse brain,Mouse liver,Mouse pancreas,Mouse thymus,Rat brain,Rat liver
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 48-177 of human IL20RB (NP_653318.2).
------	---

## 靶点信息

研究背景	
------	--

基因ID	53833
基因名	IL20RB
Swiss	Q6UXL0 ( <a href="https://www.uniprot.org/uniprotkb/Q6UXL0/entry">https://www.uniprot.org/uniprotkb/Q6UXL0/entry</a> )
别名	IL20RB,DIRS1,FNDC6,IL-20R2,IL20RB Rabbit pAb,Fibronectin type III domain containing 6

## 产品验证



Western blot analysis of IL20RB expressed in Mouse brain,Rat brain,MCF7 using IL20RB Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)