

— ABLYBIO, Help Your Research



SOX17 Rabbit pAb

货号: **AYP13319**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	PCR: Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	
实测分子量	50KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse lung,Mouse testis,Mouse kidney,Rat kidney,Rat uterus
细胞定位	nucleoplasm,nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 343-414 of human SOX17 (NP_071899.1).
------	---------------------------------------------------------------------------------------------------------------------

靶点信息

研究背景	This gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. [provided by RefSeq, Jul 2008]
基因ID	64321
基因名	SOX17
Swiss	Q9H6I2 (https://www.uniprot.org/uniprotkb/Q9H6I2/entry)
别名	VUR3,SOX17,SOX17 Rabbit pAb

产品验证



Western blot analysis of SOX17 expressed in Mouse lung, Rat kidney, 293T using SOX17 Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)