

— ABLYBIO, Help Your Research



TOB2 Rabbit pAb

货号: **AYP13341**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	27kDa/36kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,293T,MCF7,Rat skeletal muscle
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human TOB2 (N P_057356.1).
------	--

靶点信息

研究背景	TOB2 belongs to the TOB (see TOB1; MIM 605523)/BTG1 (MIM 109580) family of antiproliferative proteins, which are involved in the regulation of cell cycle progression.[supplied by OMIM, Apr 2004]
基因ID	10766
基因名	TOB2
Swiss	Q14106 (https://www.uniprot.org/uniprotkb/Q14106/entry)
别名	TOB2,APRO5,TOB4,TOBL,TROB2,TOB2 Rabbit pAb,Protein Tob4,Transducer of erbB-2 2,KIAA1663

产品验证



Western blot analysis of TOB2 expressed in Rat skeletal muscle, Jurkat, 293T using TOB2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)