

# Sall4 Rabbit pAb

货号: **AYP13352**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Mus musculus
应用	<a href="#">WB</a> <a href="#">IF/ICC</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	65kDa/112kDa
实测分子量	80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HeLa,Mouse heart
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human Sall4 (NP_065169.1).
序列	MSRRKQAKPQHINSEEDQGEQQPQQQTPEFADAAPAAPAAGELGAPVNHHPGNDEVAEDEATVKRLRREETHVCEKCC AEFFSISEFLEHKKNKCTKNPPVLMNDSEGPVPSDFSGAVLSHQPTSPGSKDCHRENGGSSSEDMKEKPD AESVVYLKTE TALPPTPQDISYLAKGKVANTNVTLQALRGTKVAVNQRSADALPAPVPGANSIPWVLEQILC

## 靶点信息

研究背景	This gene encodes a zinc finger transcription factor thought to play a role in the development of abducens motor neurons. Defects in this gene are a cause of Duane-radial ray syndrome (DRRS). Alternative splicing results in multiple transcript variants encoding different isoforms.
基因ID	57167

基因名	SALL4
Swiss	Q9UJQ4
别名	SALL4;DRRS;HSAL4;ZNF797

## 产品验证



Western blot analysis of Sall4 expressed in Jurkat, HeLa, Mouse heart using Sall4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)