

— ABLYBIO, Help Your Research



# IRAK2 Rabbit pAb

货号: **AYP13408**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Mus musculus
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	69kDa
实测分子量	75kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HeLa,HepG2,NCI-H460,Mouse liver,Mouse kidney,Rat lung
细胞定位	cytoplasm,cytosol,nucleus,plasma membrane
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human IRAK2 (NP_001561.3).
------	---

## 靶点信息

研究背景	IRAK2 encodes the interleukin-1 receptor-associated kinase 2, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. IRAK2 is reported to participate in the IL1-induced upregulation of NF-kappaB.
基因ID	3656
基因名	IRAK2
Swiss	O43187 ( <a href="https://www.uniprot.org/uniprotkb/O43187/entry">https://www.uniprot.org/uniprotkb/O43187/entry</a> )
别名	IRAK2,IRAK-2,IRAK2 Rabbit pAb

## 产品验证



Western blot analysis of IRAK2 expressed in Mouse liver,Rat lung,A-549 using IRAK2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))