

# Grp94 Rabbit pAb

货号: **AYP13456**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Homo sapiens
应用	<a href="#">WB</a> <a href="#">IHC</a> <a href="#">IF/ICC</a>
推荐浓度	<b>WB:</b> 1:500 - 1:1000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	92kDa
实测分子量	100KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,NIH/3T3,PC-12,Rat testis
细胞定位	Endoplasmic reticulum lumen,Melanosome
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 501-750 of human Grp94 (NP_003290.1).
序列	SNRTRLAKLLRFQSSHHPTDITSLDQYVERMKEKQDKIYFMAGSSRKEAESSPFVERLLKKGYEVIYLTEPVDEYCIQALPEFDGKRFQNVAKEGVKFDSEKTKESREAVEKEFEPELLNWMKDKALKDKIEKAVVSQRLTESPCALVASQYGSNGMERIMKAQAYQTGKDISTNYYASQKKTFEINPRHPLIRDMLRRIKEDDEDDKTVLIDLAVVLFETATLRSGYLLPDTKAYGDRIERMLRLSLNID

## 靶点信息

研究背景	This gene encodes a member of a family of adenosine triphosphate(ATP)-metabolizing molecular chaperones with roles in stabilizing and folding other proteins. The encoded protein is localized to melanosomes and the endoplasmic reticulum. Expression of this protein is associated with a variety of pathogenic states, including tumor formation. There is a microRNA gene located within the 5' exon of this gene. There are pseudogenes for this gene on chromosomes 1 and 15.
基因ID	7184
基因名	HSP90B1
Swiss	P14625
别名	HSP90B1;ECGP;GP96;GRP94;HEL-S-125m;HEL35;TRA1

## 产品验证



Western blot analysis of Grp94 expressed in Rat testis, MCF7, NIH/3T3 using Grp94 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)