

MAP4 Rabbit pAb

货号: **AYP13477**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	10kDa/58kDa/85kDa/91kDa/102kDa/119kDa/121kDa
实测分子量	200kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 942-1152 of human MAP4 (NP_002366.2).
序列	TENIKHQPGGGRAKVEKKTEAAATTRKPESNAVTKTAGPIASAQKQPAGKVQIVSKKVSYSHIQSKCGSKDNIKHVPGGG NVQIQNKKVDISKVSSKCGSKANIKHKPGGGDVKIESQKLNFKKAQAKVGSLDNVGHLPAGGAVKTEGGGSEAPLCPG PPAGEEPAISEAAPEAGAPTSASGLNGHPTLSGGGDQREAQTLDSQIQETSI

靶点信息

研究背景	The protein encoded by this gene is a major non-neuronal microtubule-associated protein. This protein contains a domain similar to the microtubule-binding domains of neuronal microtubule-associated protein (MAP2) and microtubule-associated protein tau (MAPT/TAU). This protein promotes microtubule assembly, and has been shown to counteract destabilization of interphase microtubule catastrophe promotion. Cyclin B was found to interact with this protein, which targets cell division cycle 2 (CDC2) kinase to microtubules. The phosphorylation of this protein affects microtubule properties and cell cycle progression. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	4134
基因名	MAP4
Swiss	P27816
别名	MAP4

产品验证



Western blot analysis of MAP4 expressed in U-87MG using MAP4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn