

OSBPL3 Rabbit pAb

货号: **AYP13577**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens,Mus musculus IF: Homo sapiens,Mus musculus
应用	WB IHC
推荐浓度	WB: 1:200 - 1:2000 IHC: 1:50 - 1:200
理论分子量	64kDa/68kDa/72kDa/93kDa/97kDa/101kDa
实测分子量	101kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HeLa,HepG2,DU145,Mouse heart,Mouse lung,Mouse brain,Rat kidney
细胞定位	cytosol,filopodium tip,nuclear membrane,perinuclear endoplasmic reticulum,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 370-520 of human OSBP L3 (NP_056365.1).
序列	SSSPSAQVIGLKNALSSALAQNITDLKERLRRHAESLLLDSPAVAKSGDNLAENSRDENRALVHQLSNESRLSITDSLSEF FDAQEVLLSPSSSENEISDDDSYVSDISDNLSDNLSDNLDNERQTLGPVLDSGREAKSRRRTCLPAPC

靶点信息

研究背景	This gene encodes a member of the oxysterol-binding protein (OSBP) family, a group of intracellular lipid receptors. Most members contain an N-terminal pleckstrin homology domain and a highly conserved C-terminal OSBP-like sterol-binding domain. The encoded protein is involved in the regulation of cell adhesion and organization of the actin cytoskeleton. Alternative splicing results in multiple transcript variants.
------	--

基因ID	26031
基因名	OSBPL3
Swiss	Q9H4L5
别名	OSBPL3;ORP-3;ORP3;OSBP3

产品验证



Western blot analysis of OSBPL3 expressed in Mouse heart,Rat kidney,Jurkat using OSBPL3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn