

— ABLYBIO, Help Your Research



# OSBPL3 Rabbit pAb

货号: **AYP13577**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Homo sapiens,Mus musculus <b>IF:</b> Homo sapiens,Mus musculus
应用	WB IHC
推荐浓度	<b>WB:</b> 1:200 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	64kDa/68kDa/72kDa/93kDa/97kDa/101kDa
实测分子量	101kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HeLa,HepG2,DU145,Mouse heart,Mouse lung,Mouse brain,Rat kidney
细胞定位	cytosol,filopodium tip,nuclear membrane,perinuclear endoplasmic reticulum,plasma membrane
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 370-520 of human OSBP L3 (NP_056365.1).
------	---

## 靶点信息

研究背景	This gene encodes a member of the oxysterol-binding protein (OSBP) family, a group of intracellular lipid receptors. Most members contain an N-terminal pleckstrin homology domain and a highly conserved C-terminal OSBP-like sterol-binding domain. The encoded protein is involved in the regulation of cell adhesion and organization of the actin cytoskeleton. Alternative splicing results in multiple transcript variants.
基因ID	26031
基因名	OSBPL3
Swiss	Q9H4L5 ( <a href="https://www.uniprot.org/uniprotkb/Q9H4L5/entry">https://www.uniprot.org/uniprotkb/Q9H4L5/entry</a> )
别名	OSBPL3,ORP-3,ORP3,OSBP3,OSBPL3 Rabbit pAb,Phosphoinositide transfer protein,KIAA0704

## 产品验证



Western blot analysis of OSBPL3 expressed in Mouse heart, Rat kidney, Jurkat using OSBPL3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))