

# RALB Rabbit pAb

货号: **AYP13597**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Mus musculus
应用	<a href="#">WB</a> <a href="#">IHC</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	23kDa/25kDa
实测分子量	23kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,HepG2,Mouse lung,Rat lung
细胞定位	Cell membrane,Cytoplasmic side,Lipid-anchor
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-206 of human RALB ( NP_002872.1).
序列	MAANKSKGQSSLALHKVIMVGGSGVGSALTLQFMYDEFVEDYEPTKADSYRKKVLDGEEVQIDILDTAGQEDYAAIRD NYFRSGEGFLLVFSITEHESFTATAEFREQILRVKAEEDKIPLLVGNKSDLEERRQVPVEEARSKAEWGVQYVETSAKTRA NVDKVFDFLMREIRTKKMSENKDKNGKKSSKNKKSFKERCCLL

## 靶点信息

研究背景	This gene encodes a GTP-binding protein that belongs to the small GTPase superfamily and Ras family of proteins. GTP-binding proteins mediate the transmembrane signaling initiated by the occupancy of certain cell surface receptors.
基因ID	5899

基因名	RALB
Swiss	P11234
别名	RALB

## 产品验证



Western blot analysis of RALB expressed in Mouse lung, Rat lung, SKOV3 using RALB Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)