

PELP1 Rabbit pAb

货号: **AYP13632**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	119kDa
实测分子量	160kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-549,HeLa,NIH/3T3,Mouse liver,Rat uterus
细胞定位	Cytoplasm,Nucleus,nucleoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1021-1180 of human PELP1 (NP_055204.3).
序列	EVEEGAPPPPTLPPALPPPEPPKVVQPEPEPEPGLLLEVEEPGTEEERGADTAPTLAPEALPSQGEVEREGESPAAGPPPQELVEEESPAPPTLLEETEDGSDKVQPPPETPAEEEMETETEAEALQEKEQDDTAAMLADFIDCPPDDEKPPPPTEPDS

靶点信息

研究背景	This gene encodes a transcription factor which coactivates transcription of estrogen receptor responsive genes and corepresses genes activated by other hormone receptors or sequence-specific transcription factors. Expression of this gene is regulated by both members of the estrogen receptor family. This gene may be involved in the progression of several types of cancer. Alternative splicing results in multiple transcript variants.
------	--

基因ID	27043
基因名	PELP1
Swiss	C9JFV4
别名	PELP1;MNAR;P160;proline

产品验证



Western blot analysis of PELP1 expressed in Mouse liver,Rat uterus,U-87MG using PELP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn