

# PPM1K Rabbit pAb

货号: **AYP13746**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Mus musculus
应用	<a href="#">WB</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	16kDa/25kDa/40kDa
实测分子量	41kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,Mouse heart,Mouse brain,Rat heart,Rat brain
细胞定位	Mitochondrion matrix
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 183-372 of human PPM1K (NP_689755.3).
序列	LLTSGTTATVALLRDGIELVVASVGDSTRAILCRKGGKPMKLTIDHTPERKDEKERIKKCGGFVAWNSLQPHVNGRLAMTRSI GDLDLKTSGVIAEPETKRIKLHHADDSFLVLTDDGINFMVNSQEICDFVNQCHDPNEAAHAVTEQAIQYGTEDNSTAVVVP FGAWGKYKNSEINFSSRSFASSGRWA

## 靶点信息

研究背景	This gene encodes a member of the PPM family of Mn <sup>2+</sup> /Mg <sup>2+</sup> -dependent protein phosphatases. The encoded protein, essential for cell survival and development, is targeted to the mitochondria where it plays a key role in regulation of the mitochondrial permeability transition pore.
基因ID	152926

基因名	PPM1K
Swiss	Q8N3J5
别名	PPM1K;BDP;MSUDMV;PP2Ckappa;PP2Cm;PTMP;UG0882E07

## 产品验证



Western blot analysis of PPM1K expressed in Mouse heart,Rat heart,A-431 using PPM1K Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane . Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)