

— ABLYBIO, Help Your Research



# ZNF76 Rabbit pAb

货号: **AYP13752**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Homo sapiens
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	56kDa/61kDa
实测分子量	61kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse testis,Rat serum
细胞定位	nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-170 of human ZNF76 (NP_003418.2).
------	---

## 靶点信息

研究背景	May be involved in transcriptional regulation.
------	--

基因ID	7629
基因名	ZNF76
Swiss	P36508 ( <a href="https://www.uniprot.org/uniprotkb/P36508/entry">https://www.uniprot.org/uniprotkb/P36508/entry</a> )
别名	ZNF76,D6S229E,ZNF523,Zfp523,ZNF76 Rabbit pAb,Zinc finger protein 523

## 产品验证



Western blot analysis of ZNF76 expressed in HeLa, Mouse testis, Rat serum using ZNF76 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)