

# CSF2RA Rabbit pAb

货号: **AYP13761**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Mus musculus
应用	<a href="#">WB</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	26kDa/31kDa/32kDa/38kDa/43kDa/46kDa/49kDa
实测分子量	80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,DU145,HepG2,Mouse spleen,Mouse lung
细胞定位	Cell membrane,Secreted,Single-pass type I membrane protein
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 23-320 of human CSF2R A (NP_758448.1).
序列	EKSDLRTVAPASSLNVRFDSTMTNLSWDCQENTTFSKCFLTDKKNRVVEPRLSNNECSCTFREICLHEGVTFEVHVNTSQ RGFQQKLLYPNSGREGTAAQNFSCFIYNADLMNCTWARGPTAPRDVQYFLYIRNSKRREIRCPYYIQDSGTHVGCHLDN LSGLTSRNYFLVNGTSREIGIQFFDSSLDTKKIERFNPPSNVTVRCNTTHCLVRWKQPRTYQKLSYLDFOYQLDVHRKNTQ PGTENLLINVSGDLENRYNFPSSSEPRAKHSVKIRAADVRLNWSWSEAIEFGSDDG

## 靶点信息

研究背景	The protein encoded by this gene is the alpha subunit of the heterodimeric receptor for colony stimulating factor 2, a cytokine which controls the production, differentiation, and function of granulocytes and macrophages. The encoded protein is a member of the cytokine family of receptors. This gene is found in the pseudoautosomal region (PAR) of the X and Y chromosomes. Multiple transcript variants encoding different isoforms have been found for this gene, with some of the isoforms being membrane-bound and others being soluble.
基因ID	1438
基因名	CSF2RA
Swiss	P15509
别名	CSF2RA;CD116;CDw116;CSF2R;CSF2RAX;CSF2RAY;CSF2RX;CSF2RY;GM-CSF-R-alpha;GMCSFR;GMR;SMDP4

## 产品验证



Western blot analysis of CSF2RA expressed in Mouse spleen, HeLa, DU145 using CSF2RA Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)