

# CDIP1 Rabbit pAb

货号: **AYP13931**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Homo sapiens
应用	<a href="#">WB</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	13kDa/17kDa/21kDa
实测分子量	22kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Rat heart
细胞定位	nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CDIP1 (NP_001185983.1).
序列	YTPGPYPGPGGHTATVLPVSGAATTVTVLQGEIFEGAPVQVCPHCQQAITTKISYEIGLMNFVLGFFCCFMGCDLGCCLIPCLINDFKDVTHTCPSCKAY

## 靶点信息

研究背景	Acts as an important p53/TP53-apoptotic effector. Regulates TNF-alpha-mediated apoptosis in a p53/TP53-dependent manner.
基因ID	29965
基因名	CDIP1

Swiss	Q9H305
别名	CDIP1;C16orf5;CDIP;I1;LITAFL

## 产品验证



Western blot analysis of CDIP1 expressed in Mouse heart,Rat heart using CDIP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)