

— ABLYBIO, Help Your Research



CDIP1 Rabbit pAb

货号: **AYP13931**

产品信息

| | |
|-------|--|
| 反应 | Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | WB: Homo sapiens |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:2000 |
| 理论分子量 | 13kDa/17kDa/21kDa |
| 实测分子量 | 22kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Mouse heart,Rat heart |
| 细胞定位 | nucleus |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|---|
| 抗原信息 | A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CDIP1 (NP_001185983.1). |
|------|---|

靶点信息

| | |
|------|--|
| 研究背景 | Acts as an important p53/TP53-apoptotic effector. Regulates TNF-alpha-mediated apoptosis in a p53/TP53-dependent manner. |
|------|--|

| | |
|-------|---|
| 基因ID | 29965 |
| 基因名 | CDIP1 |
| Swiss | Q9H305 (https://www.uniprot.org/uniprotkb/Q9H305/entry) |
| 别名 | CDIP1,C16orf5,CDIP,I1,LITAF,CDIP1 Rabbit pAb,Cell death involved p53-target,Cell death-inducing protein,LITAF-like protein,Lipopolysaccharide-induced tumor necrosis factor-alpha-like protein,Transmembrane protein I1 |

产品验证



Western blot analysis of CDIP1 expressed in Mouse heart, Rat heart using CDIP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)