

— ABLYBIO, Help Your Research



SERBP1 Rabbit pAb

货号: **AYP13932**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	42kDa/43kDa/44kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,A-549,HeLa,Mouse lung,Rat liver,Rat lung
细胞定位	Cytoplasm,Nucleus,perinuclear region
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-114 of human SERBP1 (NP_001018078.1).
------	---

靶点信息

研究背景	May play a role in the regulation of mRNA stability. Binds to the 3'-most 134 nt of the SERPINE1/PAI1 mRNA, a region which confers cyclic nucleotide regulation of message decay. Seems to play a role in PML-nuclear bodies formation.
基因ID	26135
基因名	SERBP1
Swiss	Q8NC51 (https://www.uniprot.org/uniprotkb/Q8NC51/entry)
别名	SERBP1, CGI-55, CHD3IP, HABP4L, PAI-RBP1, PAIRBP1, SERBP1 Rabbit pAb, PAI1 RNA-binding protein 1, Plasminogen activator inhibitor 1 RNA-binding protein

产品验证



Western blot analysis of SERBP1 expressed in Mouse lung, Rat liver, U-251MG using SERBP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)