

# SLC25A17 Rabbit pAb

货号: **AYP13935**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Homo sapiens
应用	<a href="#">WB</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	34kDa
实测分子量	34kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Mouse kidney
细胞定位	Cytoplasm,Multi-pass membrane protein,Peroxisome membrane
纯化	Affinity purification

## 抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human SLC25A17 (NP_06349.1).
序列	IDAFHQIIRDEGISALWNGTFPSSLVFNPAIQMFYEGLKRQLLKKRMKLSLDVFIIGAVAKAIATTVTYPLQTVQSILRFGR HRLNPENRTLGLSLRNI

## 靶点信息

研究背景	This gene encodes a peroxisomal membrane protein that belongs to the family of mitochondrial solute carriers. It is expressed in the liver, and is likely involved in transport. Alternative splicing results in multiple transcript variants.
基因ID	10478

基因名	SLC25A17
Swiss	O43808
别名	SLC25A17;PMP34

## 产品验证



Western blot analysis of SLC25A17 expressed in Mouse liver, Mouse kidney using SLC25A17 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)