

# WNT3A Rabbit pAb

货号: **AYP14019**

## 产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>IHC:</b> Homo sapiens
应用	<a href="#">IHC</a>
推荐浓度	<b>IHC:</b> 1:50 - 1:200
理论分子量	39kDa/42kDa
实测分子量	39kDa/42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-352 of human WNT3A (NP_149122.1).
序列	MAPLGYFLLLC SLKQALGSYPIWWSLAVGPQYSSLGSQPILCASIPGLVPKQLRFCRNYVEIMPSVAEGIKIGIQECQHQFRG RRWNCTTVHDSLAIFGPVLDKATRESAFVHAIASAGVAFVTRSCAEGTAAICGCSRHQSGSPGKGWKWGGCSEIEFG GMVSREFADARENRPDARSAMNRHNN EAGRQAIASHMHLKCKCHGLSGSCEVKTCWWSQPDFRAIGDFLKDKYDSAS EMVVEKHRESRGWVETLRPRYTYFKVPTERDLVYEEASP NFCEPNPETGSGFTRDRTCNVSSHGIDGCDLLCCGRGHNAR AERRREKCRCVFHWCCYVSCQECTRVYDVHTCK

## 靶点信息

研究背景	The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It encodes a protein which shows 96% amino acid identity to mouse Wnt3A protein, and 84% to human WNT3 protein, another WNT gene product. This gene is clustered with WNT14 gene, another family member, in chromosome 1q42 region.
基因ID	89780
基因名	WNT3A
Swiss	P56704
别名	WNT3A

#### 产品验证

#### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)