

— ABLYBIO, Help Your Research



PELP1 Rabbit pAb

货号: **AYP14033**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens
应用	IHC
推荐浓度	IHC: 1:50 - 1:200
理论分子量	119kDa
实测分子量	119kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm,Nucleus,nucleoplasm
纯化	Affinity purification

抗原信息

抗原信息	Recombinant protein of human PELP1
------	------------------------------------

靶点信息

研究背景	This gene encodes a transcription factor which coactivates transcription of estrogen receptor responsive genes and corepresses genes activated by other hormone receptors or sequence-specific transcription factors. Expression of this gene is regulated by both members of the estrogen receptor family. This gene may be involved in the progression of several types of cancer. Alternative splicing results in multiple transcript variants.
基因ID	27043
基因名	PELP1
Swiss	Q8IZL8 (https://www.uniprot.org/uniprotkb/Q8IZL8/entry)
别名	PELP1,MNAR,P160,proline,PELP1 Rabbit pAb,Modulator of non-genomic activity of estrogen receptor,Transcription factor HMX3,HMX3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)