

USP30 Rabbit pAb

货号: **AYP14079**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Homo sapiens
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	58kDa
实测分子量	71kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Rat lung
细胞定位	Mitochondrion outer membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 57-390 of human USP30 (NP_116052.2).
序列	TERKKRRKGLVPLVNLGNTCFMNSLLQGLSACPAFIRWLEEFSTQYSRDQKEPPSHQYLSLTLHLLKALSCQEVTDDEV LDASCLLDVLRMYRWQISSFEEQDAHELPHVITSSLEDERDRQPRVTHLFDVHSLEQQSEITPKQITCRTRGSPHPTS KSHQPFHGRLTSMVCKHCEHQSPVRFDTFDSL SLSIPAATWGHPLTLDHCLHHFISSESVRDVVCDNCTKIEAKGTLNG EKVEHQRTTFVKQLKGLPQCLCIHLQRLSWSSHGTPPKRHEHVQFNEFLMMDIYKYHLLGHKPSQHNPKNKPNPGPTL ELQDGPAPTPTV

靶点信息

研究背景	USP30, a member of the ubiquitin-specific protease family (see USP1, MIM 603478), is a novel mitochondrial deubiquitinating (DUB) enzyme (Nakamura and Hirose, 2008 [PubMed 18287522]).[supplied by OMIM, Dec 2008]
------	---

基因ID	84749
基因名	USP30
Swiss	Q70CQ3
别名	USP30

产品验证



Western blot analysis of USP30 expressed in Mouse lung,Rat lung using USP30 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn