

# NOXO1 Rabbit pAb

货号: **AYP14217**

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Mus musculus
应用	<a href="#">WB</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	40kDa/41kDa
实测分子量	41kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Mouse small intestine, Mouse stomach
细胞定位	Cell membrane, Cytoplasmic side, Peripheral membrane protein
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 121-370 of human NOX O1 (NP_653204.1).
序列	QPLDLEPALPPGSRVILPTPEEQPLSRAAGRLSIHSLEAQLRCLQPFCTQDTRDRPFQAQAQESLDVLLRHPSGWWLVEN EDRQTAWFPAPYLEEAAPGQGREGGPSLGSSGPOFCASRAYESSRADELSPAGARVRVLETSDRGWWLCRYGDRAGL LPAVLLRPEGLGALLSGTGFRGGDDPAGEARGFPEPSQATAPPPTVPTRPSPGAIQSRCCTVTRRALERRPRRQGRPRGCV DSVPHPTTEQ

## 靶点信息

研究背景	This gene encodes an NADPH oxidase (NOX) organizer, which positively regulates NOX1 and NOX3. The protein contains a PX domain and two SH3 domains. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	124056

基因名	NOXO1
Swiss	Q8NFA2
别名	NOXO1;P41NOX;P41NOXA;P41NOXB;P41NOXC;SH3PXD5;SNX28

## 产品验证



Western blot analysis of NOXO1 expressed in Mouse small intestine, Mouse stomach using NOXO1 Rabbit p Ab at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)