

— ABLYBIO, Help Your Research



# NOXO1 Rabbit pAb

货号: **AYP14217**

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Mus musculus
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	40kDa/41kDa
实测分子量	41kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Mouse small intestine, Mouse stomach
细胞定位	Cell membrane, Cytoplasmic side, Peripheral membrane protein
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 121-370 of human NOXO1 (NP_653204.1).
------	---

## 靶点信息

研究背景	This gene encodes an NADPH oxidase (NOX) organizer, which positively regulates NOX1 and NOX3. The protein contains a PX domain and two SH3 domains. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	124056
基因名	NOXO1
Swiss	Q8NFA2 ( <a href="https://www.uniprot.org/uniprotkb/Q8NFA2/entry">https://www.uniprot.org/uniprotkb/Q8NFA2/entry</a> )
别名	NOXO1,P41NOX,P41NOXA,P41NOXB,P41NOXC,SH3PXD5,SNX28,NOXO1 Rabbit pAb,NADPH oxidase regulatory protein,Nox organizer 1,Nox-organizing protein 1,SH3 and PX domain-containing protein 5

## 产品验证



Western blot analysis of NOXO1 expressed in Mouse small intestine, Mouse stomach using NOXO1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)