

# FSTL3 Rabbit pAb

货号: AYP14230

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Homo sapiens
应用	<a href="#">WB</a>
推荐浓度	<b>WB:</b> 1:1000 - 1:2000
理论分子量	24kDa/27kDa
实测分子量	27kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,DU145,MCF7,SKOV3,Mouse lung,Mouse kidney,Rat lung
细胞定位	Nucleus,Secreted
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 27-263 of human FSTL3 (NP_005851.1).
序列	MGSGNPAPGGVCWLQQGQEATCSLVLQTDVTRAECASGNIDTAWSNLTHPGNKINLLGFLGLVHCLPCKDSCDGVEC GPGKACRMLGGRPRCECAPDCSGLPARLQVCGSDGATYRDECELRAARCRGHPDLSVMYRGRCRKSCEHVVCPRPQSC VVDQTGSAHCVVCRAPCPVPSSPGQELCGNNAVITYISSCHMRQATCFLGRSIGVRHAGSCAGTPEEPPGGESAEEEENF V

## 靶点信息

研究背景	Follistatin-like 3 is a secreted glycoprotein of the follistatin-module-protein family. It may have a role in le ukemogenesis.
基因ID	10272

基因名	FSTL3
Swiss	O95633
别名	FSTL3;FLRG;FSRP

## 产品验证



Western blot analysis of FSTL3 expressed in Mouse lung,Rat lung,HeLa using FSTL3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)