

— ABLYBIO, Help Your Research



RETNLB Rabbit pAb

货号: **AYP14390**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	10kD,25-40kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-111 of human RETNLB (NP_115968.1).
------	--

靶点信息

研究背景	Probable hormone.
------	-------------------

基因ID	84666
基因名	RETNLB
Swiss	Q9BQ08
别名	RETNLB,FIZZ1,FIZZ2,HXCP2,RELM-beta,RELMb,RELMbeta,XCP2,RETNLB Rabbit pAb,Colon and small intestine-specific cysteine-rich protein,Colon carcinoma-related gene protein,Cysteine-rich secreted protein A1 2-alpha-like 1,Cysteine-rich secreted protein FIZZ2,CCRG,RETNL2

产品验证



Western blot analysis of RETNLB expressed in HUVEC, mouse heart using RETNLB Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)