

— ABLYBIO, Help Your Research



MIB2 Rabbit pAb

货号: **AYP14436**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	Co-IP: Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,BxPC-3,U-87MG,Mouse lung,Mouse testis,Mouse brain,Rat brain,Rat testis
细胞定位	cytoplasm,cytosol,endosome
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 300-580 of human MIB2 (NP_001164158.1).
------	---

靶点信息

研究背景	The protein encoded by this gene is an E3 ubiquitin protein ligase that mediates ubiquitination of proteins in the Notch signaling pathway. The encoded protein may be a suppressor of melanoma invasion. [provided by RefSeq, Mar 2017]
基因ID	142678
基因名	MIB2
Swiss	Q96AX9 (https://www.uniprot.org/uniprotkb/Q96AX9/entry)
别名	ZZANK1,ZZZ5,MIB2,MIB2 Rabbit pAb,Mind bomb homolog 2,Novel zinc finger protein,Putative NF-kappa-B-activating protein 002N,RING-type E3 ubiquitin transferase MIB2,Skeletrophin,Zinc finger ZZ type with a nkyrin repeat domain protein 1,SKD

产品验证



Western blot analysis of MIB2 expressed in Mouse lung,Rat brain,HeLa using MIB2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)