

# MIB2 Rabbit pAb

货号: **AYP14436**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>Co-IP:</b> Mus musculus
应用	<a href="#">WB</a> <a href="#">IF/ICC</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,BxPC-3,U-87MG,Mouse lung,Mouse testis,Mouse brain,Rat brain,Rat testis
细胞定位	cytoplasm,cytosol,endosome
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 300-580 of human MIB2 (NP_001164158.1).
序列	ASVTWADGTTNVYRVGHKGKVDLKCVGGEAAGGFYKDHLPRLGKPAELQRRVSADSQPFQHGDKVKCLLDTDLVREM QEGHGGWNPRMAEFIGQTGTVHRITDRGDVRVQFNHETRWFHPGALTKHHSFWVGDVVRVIGDLDTVKRLQAGHGE WTDDMAPALGRVGVKVFVDGNLRVAVAGQRWTFSPSCLVAYRPEEDANLDAERARENKSSLSVALDKLRAQKSDP EHPGRLVVEVALGNAARALDLLRRRPEQVDTKNQGRALQVAAYLGQVE

## 靶点信息

研究背景	The protein encoded by this gene is an E3 ubiquitin protein ligase that mediates ubiquitination of proteins in the Notch signaling pathway. The encoded protein may be a suppressor of melanoma invasion. [provided by RefSeq, Mar 2017]
------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

基因ID	142678
基因名	MIB2
Swiss	Q96AX9
别名	ZZANK1;ZZZ5;MIB2

## 产品验证



Western blot analysis of MIB2 expressed in Mouse lung,Rat brain,HeLa using MIB2 Rabbit pAb at 1:1000 . Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)