

— ABLYBIO, Help Your Research



# PML/RARA Rabbit pAb

货号: **AYP14921**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>ChIP:</b> unknow
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	47-48kDa/62-97kDa/39kDa/50kDa
实测分子量	51kDa/130kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,NIH/3T3,C6
细胞定位	cytoplasm,cytosol,nuclear matrix,nuclear membrane,nucleolus,nucleoplasm,nucleus,PML body
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant protein of human PML/RARA
------	---------------------------------------

## 靶点信息

研究背景	
基因ID	5371,5914

基因名	PML,RARA
Swiss	P29590 ( <a href="https://www.uniprot.org/uniprotkb/P29590/entry">https://www.uniprot.org/uniprotkb/P29590/entry</a> ), P10276 ( <a href="https://www.uniprot.org/uniprotkb/P10276/entry">https://www.uniprot.org/uniprotkb/P10276/entry</a> )
别名	PML/RARA,PML/RARA Rabbit pAb,PML,RARA,E3 SUMO-protein ligase PML,Promyelocytic leukemia protein,RING finger protein 71,RING-type E3 SUMO transferase PML,Tripartite motif-containing protein 19,Nuclear receptor subfamily 1 group B member 1,MYL,PP8675,RNF71,TRIM19

## 产品验证



Western blot analysis of PML/RARA expressed in 293T,NIH/3T3,C6 using PML/RARA Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)