

DTX2 Rabbit pAb

货号: **AYP14922**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: 293T
应用	WB IHC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200
理论分子量	62kDa/67kDa
实测分子量	67KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human DTX2 (NP_065943.2).
序列	MAMAPSPSLVQVYTSPAAVAVWEWQDGLGTWHPYSATVCSFIEQQFVQQKGQRFGLGSLAHSIPLGQADPSLAPYIIDLPSWTQFRQDTGTMRAVRRHLFPQHSAPGRGVVWEWLSDDGSWTAYEASVCDYLEQQVARGNQLVDLAPLGYNITVNYTTHTQTNKTSSFCRSVRRQAGPPYPVTIIAPPGHTGVACSCHQCLSGSRTGPVSGRGRHSMTNLPAYPVPQHPPHRTASVFGTHQAFAPYINKPSLSGARSAPRLNTTNAWGAAPPSLGSQPLY

靶点信息

研究背景	DTX2 functions as an E3 ubiquitin ligase (Takeyama et al., 2003 [PubMed 12670957]).[supplied by OMIM, Nov 2009]
基因ID	113878

基因名	DTX2
Swiss	Q86UW9
别名	DTX2;RNF58

产品验证



Western blot analysis of DTX2 expressed in A-549 using DTX2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn