

— ABLYBIO, Help Your Research



# FAH Rabbit pAb

货号: **AYP14927**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> HepG2,H460
应用	WB IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:10 - 1:100
理论分子量	38kDa/46kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NCI-H460,HepG2,MCF7,Raji
细胞定位	cytosol,extracellular exosome
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-419 of human FAH (NP_000128.1).
------	---

## 靶点信息

研究背景	This gene encodes the last enzyme in the tyrosine catabolism pathway. FAH deficiency is associated with Type 1 hereditary tyrosinemia (HT).
基因ID	2184
基因名	FAH
Swiss	P16930 ( <a href="https://www.uniprot.org/uniprotkb/P16930/entry">https://www.uniprot.org/uniprotkb/P16930/entry</a> )
别名	FAH,FAH Rabbit pAb,Beta-diketonase,Fumarylacetoacetate hydrolase

## 产品验证



Western blot analysis of FAH expressed in NCI-H460,HepG2,MCF7 using FAH Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))