

FAM160B2 Rabbit pAb

货号: **AYP15061**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:500 - 1:1000
理论分子量	82kDa
实测分子量	82kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,A-549,SH-SY5Y,Mouse testis,Rat testis
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 74-303 of human FAM160B2 (NP_073586.5).
序列	GEAGPCLEYLLQHKILETLCTLGKAEYPPGMRQQVFQFFSKVLAQVQHPLLHYLSVHRPVQKLLRLGGTASGSVTEKEEVQFTTVLCSKIQDPELLAYILEGKKIVGRKKACGEPTALPKDTTSHGDKDCSHDGAPARPQLDGESCGAQLNSHMPAETEE LDGGTTESNLITSLGLCQSKSRVALKAQENLLLLVSMASPAATYLVQSSACCPAIVRHLCQLYRS

靶点信息

研究背景	Able to activate MAPK/ERK and TGFB signaling pathways. May regulate the activity of genes involved in intestinal barrier function and immunoprotective inflammation (By similarity. May play a role in cell proliferation).
------	---

基因ID	64760
基因名	FAM160B2
Swiss	Q86V87
别名	FAM160B2; RAI16; protein FAM160B2

产品验证



Western blot analysis of FAM160B2 expressed in Mouse testis, Rat testis, 293T using FAM160B2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30u g per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn