

— ABLYBIO, Help Your Research



Phospho- α -Synuclein-S129 Rabbit pAb

货号: **AYP15067**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	11kDa/13kDa/14kDa
实测分子量	20kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Cell junction,Cytoplasm,Membrane,Nucleus,Secreted,cytosol,synapse
纯化	Affinity purification


抗原信息

抗原信息	A synthetic phosphorylated peptide around S129 of human alpha-Synuclein (NP_000336.1).
------	--

靶点信息

研究背景	Alpha-synuclein is a member of the synuclein family, which also includes beta- and gamma-synuclein. Synucleins are abundantly expressed in the brain and alpha- and beta-synuclein inhibit phospholipase D2 selectively. SNCA may serve to integrate presynaptic signaling and membrane trafficking. Defects in SNCA have been implicated in the pathogenesis of Parkinson disease. SNCA peptides are a major component of amyloid plaques in the brains of patients with Alzheimer's disease. Alternatively spliced transcripts encoding different isoforms have been identified for this gene.
基因ID	6622
基因名	SNCA
Swiss	P37840 (https://www.uniprot.org/uniprotkb/P37840/entry)
别名	PD1,NACP,PARK1,PARK4,SNCA,Phospho- α -Synuclein-S129 Rabbit pAb,Non-A beta component of AD amyloid,Non-A4 component of amyloid precursor

产品验证



Western blot analysis of Phospho- α -Synuclein-S129 expressed in Mouse brain using Phospho- α -Synuclein-S129 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)