

# ZFP42 Rabbit pAb

货号: **AYP15119**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	34kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ZFP42 (NP_777560.2).
序列	MSQQLKKRAKTRHQKGLGGRAPSGAKPRQKGSSQDLQAEIEPVSAVWALCDGYVCYEPGPQALGGDDFSDCYIECVIR GEFSQPILEEDSLFESLEYLKKGSEQLSQKVFEASSLECSLEYMKKGVKKELPQKIVGENS

## 靶点信息

研究背景	Involved in the reprogramming of X-chromosome inactivation during the acquisition of pluripotency. Required for efficient elongation of TSIX, a non-coding RNA antisense to XIST. Binds DXPas34 enhancer within the TSIX promoter. Involved in ES cell self-renewal (By similarity).
基因ID	132625

基因名	ZFP42
Swiss	Q96MM3
别名	ZFP42;REX-1;REX1;ZNF754;zfp-42

## 产品验证



Western blot analysis of ZFP42 expressed in 293T using ZFP42 Rabbit pAb at 1:1000. Secondary antibody : HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)